

Performance Report on Future Fields' Recombinant Human FGF2

Empowering Life Science Innovators

At Future Fields, we are passionate about helping you achieve your research goals. We know that you are constantly facing new challenges, and unreliable recombinant proteins that break the bank should not be one of them. That is why we have designed our products with quality, affordability, and sustainability in mind—to pass on these benefits to you, the scientist.

In this piece, we lay out an in-depth performance analysis of our latest product, [Recombinant Human FGF2](#), and demonstrate its ability to stimulate proliferation in a dose-dependent manner, all the while keeping the cell morphology and other characteristics intact.

Here is the short version:

1. **Performance:** Functionality of Human FGF2 validated using multiple metrics
 - 1.1. SPR readout confirming protein identity and binding
 - Ligand-Receptor interaction (< 5 nM binding affinity) and
 - Ligand-Antibody interaction (sub-nanomolar binding affinity)
 - 1.2. Consistent EC₅₀ ~ 1-3 ng/ml in cell-based proliferation assay
 - 1.3. Fluorescent images to demonstrate healthy cell morphology
2. **Purity:** >95 % purity
3. **Safety:**
 - 3.1. Negative for Mycoplasma contamination
 - 3.2. Below safety threshold (< 1 EU/ug) for Endotoxin

Let us dive in.

1 Performance

First, let us start at the heart of our product—its performance. We use multiple methods to qualify our products; biophysical techniques like Surface Plasmon Resonance (SPR) can tell us how tight the binding of the ligand to the substrate (receptor, ligand-specific antibody or other known interactors), as well as cell-based functionality studies to determine that our products work as expected, and they meet your standards for cell culture applications.

1.1 Surface Plasmon Resonance (SPR)

Using FGF-receptor (human), or anti-FGF2 antibody immobilized on a Carboxymethyl Dextran chip, we demonstrate Future Fields' Human FGF2 has an affinity constant of <1 nM on antibody, and ~ 5 nM on the receptor (Fig. 1). This **nanomolar affinity** indicates a very tight ligand-substrate interaction.

Fig. 1 - Graph of SPR Sensorgram




Figure 1. SPR sensorgram showing tight affinity of Human FGF2 on FGF-receptor bound to the substrate. Curve-fit indicates $K_d < 5$ nM for ligand-receptor interaction.

1.2 Cell-Based Studies (Various Assays)

1.2.1 Proliferation

To characterize Future Fields' Human FGF2 for their functionality, we use multiple cell-based assays. One among them is to assess FGF2 ability to drive proliferation in NIH-3T3 cells (mouse-derived, to demonstrate cross-reactivity). Using a metabolic activity-based readout, we determine that FGF2 is able to drive proliferation in the cells in a dose-dependent manner ($EC_{50} = 1-3$ ng/ml), and is highly comparable to existing competitors' offerings. We also demonstrate lot-to-lot consistency of proteins via cell proliferation, to reduce variability concerns for your research purposes (Fig. 2).



Fig. 2 - % of Proliferation of NIH-3T3 Cells Stimulated with Human FGF2




Figure 2. Graph represents multiple lots of Future Fields' Human FGF2 tested for their ability to drive proliferation of NIH-3T3 cells, using alamar blue assay, in comparison to competitor standard. EC₅₀ is determined to be in the low nanomolar range indicating active ligands.

1.3 Visuals Demonstrating Cell Morphology

1.3.1 Brightfield Cell Images and Population Doubling Time

These representative images (Fig. 3) show cell number increase in time (population doubling time ~20 hours), as well as cell morphology including shape, size, and intact nuclei that affirm that cells at low and high density are **healthy and proliferating robustly**.





Fig. 3 - Graph Showing Cell Proliferation and Images Tracked in Time




Figure 3. Graph representing NIH-3T3 cells proliferation in time, when cultured with Human FGF2. Population doubling time was determined to be ~20 hours, which is in line with literature.

1.3.2 Fluorescent Cell Images

In addition to the above images, we also highlight intact nuclei (DAPI stained) and actin filament structures (phalloidin stained) to demonstrate the cell boundaries and expected morphology of healthy cells that are growing in media rich with Future Fields' Human FGF2 (Fig. 4).

Fig. 4 - Fluorescent Cell Images




Figure 4. (1) Representative merged images of actin and phalloidin stained NIH-3T3 cells in active state of proliferation and growth. (2) DAPI staining indicates intact nuclei structures. (3) Phalloidin staining indicates actin filaments, showing cell boundaries and demarcations.





2 Purity

We purify our growth factors in a proprietary column-based purification strategy (Fig. 5) that helps maintain the structure and functionality of the product throughout the process, while also ensuring that no or limited host components are part of the final products. In Fig. 5, we show **purity > 95 %** of Future Fields' Human FGF2 from different lots in comparison to industry standard.

Fig. 5 - SDS-PAGE Gel Showing Purity of Product




Figure 5. Representative image of SDS-PAGE showing different lots of Future Fields' FGF2 in comparison to competitor standard; proteins run at expected size and no contaminant proteins are visible indicating high purity proteins.

3 Safety

We understand safety is paramount when it comes to cell culture studies. Therefore, we performed assays for the two common contaminants in cell culture media.

3.1 *Mycoplasma*

A primary concern is *Mycoplasma*, the intracellular antibiotic-resistant contaminant that can have numerous effects on mammalian cell culture. Reassuringly, our proteins are always tested for several species of *Mycoplasma*, *Acholeplasma*, and *Ureaplasma* using a sensitive PCR-based test (Fig. 6) that can detect 1–5 fg of *Mycoplasma* DNA corresponding to 2–5 *Mycoplasma* per sample volume; this ensures we send products that are **negative** for these contaminations.

Fig. 6 - PCR Gel of Mycoplasma Test Results




Figure 6. Representative PCR gel image showing positive control for mycoplasma band at 250 bp, and internal control at 191 bp, and negative mycoplasma bands for all lots tested indicating no mycoplasma contamination.



3.2 Endotoxin

Our approach to quality includes rigorous testing for endotoxin, a common name for Lipopolysaccharides (LPS), derived from bacteria that can trigger apoptotic pathways in mammalian cell cultures. Future Fields' Human FGF2 has consistently shown (Fig. 7) to be **below the threshold of FBS (<1 EU/ug)**, thereby reducing any fears of unwanted deleterious cell signaling and effects.

Fig. 7 - Lot-to-Lot Endotoxin Test Results




Figure 7. Endotoxin was measured using Chromogenic assay, and all lots tested showed below established threshold values for endotoxin levels indicating safety for mammalian cell cultures.

In conclusion, the above characteristics of Future Fields' Recombinant Human FGF2 produced by the EntoEngine™ can meet your expectations for consistency of high performance and reliability of supply chain to support your research and production goals.



Designing Proteins with Quality, Affordability, and Sustainability in Mind

We're a [green-certified](#), bioreactor-free lab. Our fly-based [EntoEngine™ platform](#) produces complex recombinant proteins more sustainably and 30 times faster than legacy systems in steel tanks. That's why we can sell our products at a fraction of the cost of industry prices today. You can read more about our story in [Forbes](#), penned by the CEO of SynBioBeta.

Our recombinant proteins are also incredibly versatile. Whether you're studying disease pathways, developing therapeutics, or exploring the frontiers of synthetic biology, our recombinant proteins are the reliable, high-performance solution you can count on. We're constantly expanding our product offerings to meet the evolving needs of the biotech industry, and [we're always happy to work with you](#) to meet your specific needs.

With Future Fields, you're not just buying a product—you're partnering with [a community of forward-thinking researchers](#) who are committed to pushing the limits of what's possible in biotech. We work for humanity and the planet, which is why [1 % of your purchase proceeds](#) will go towards environmental initiatives around the world.

We're not just a vendor—we're a partner in your research success and your champion for positive change.

Get started by [requesting a sample today](#).

